

When printing, make sure that "Page Scaling" is set to "none".
To verify that the file has been printed correctly, this line should measure 8 1/2".

GeeWhiz® External Condom Catheter Sizing Guides

Order Online at:
www.seekwellness.com/mystore
For more information or to order by phone, call (800) 840-9301 or email customer.service@seekwellness.com


Cut out Sizing Strip and wrap **SNUGGLY** around shaft of penis to find correct condom size. If arrow tip touches black square between sizes, use next smaller size.

Condom Size in MM

Condom Size (MM)	Condom Size (Inches)
29	3.59"
32	3.96"
36	4.45"
41	5.07"

Use Rounded Sizing Guides or Sizing Guide Strip to Select Correct Condom Size
Estimate correct size by placing rounded guides **over shaft** of penis behind the head of penis to check size.

For best results select size that fits SNUGGLY over the shaft of the penis.



PR-L0070-041908